

Palladin-Specific antibody



Catalog Number: 113563

Product name

Palladin-Specific antibody

Specificity

Human; other species not tested.

Antibody description

Palladin-Specific Rabbit Polyclonal antibody. Positive WB detected in HeLa cells. Positive IP detected in HeLa cells. Observed molecular weight by Western-blot: 120kd

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_016081). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:1000

Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:113563(PALLD,palladin-Specific antibody) at dilution of 1:200



IP Result of anti-PALLD,palladin-Specific (IP:Catalog No:113563, 4ug; Detection:Catalog No:113563 1:200) with HeLa cells lysate 1500ug.