

Anti-GSK3 beta (2E6) Mouse antibody



Catalog Number: 168263

Product name

Anti-GSK3 beta (2E6) Mouse antibody

Specificity

Human, Mouse, Rat

Antibody description

GSK3 beta (2E6) Mouse monoclonal antibody

Preparation

Antigen: Purified recombinant GSK-3 β protein fragments expressed in E.coli.

Formulation

PBS(pH 7.4) containing with 0.03% Proclin300 and 50% glycerol.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Mouse IgG2a

Applications


WB, ICC/IF, IP, FC

Dilutions


WB: 1/1000

;ICC: 1/200


Validations



Western blot analysis of extracts from A549, HeLa, 3T3 and Rat Brain cell lysates using GSK-3 beta mouse mAb(1:1000 diluted). Predicted band size:46KDa. Observed band size:46KDa.



Immunocytochemistry staining of HeLa cells fixed with -20°C Methanol and using GSK-3 beta mouse



mAb (dilution 1:200).

Immunoprecipitation analysis of HeLa cell lysates using GSK3 beta (2E6) Mouse mAb. GSK-3 beta Rabbit mAb was used for the western blot analysis (1:1000 diluted).