

# CEP110 antibody



Catalog Number: 109263

## Product name

CEP110 antibody

## Immunogen

[Human CEP110 Recombinant protein \(His tag\)](#)

## Specificity

Human; other species not tested.

## Antibody description

CEP110 Rabbit Polyclonal antibody. Positive WB detected in Transfected HEK-293 cells.

## Preparation

This antibody was obtained by immunization of CEP110 recombinant protein (Accession Number: NM\_007018). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:2000

## Validations



Transfected HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:109263(CEP110 Antibody) at dilution of 1:700