

Cyclin B2 antibody



Catalog Number: 109663

Product name

Cyclin B2 antibody

WB: 1:500-1:5000

IP: 1:500-1:5000

Specificity

Human, Mouse, Rat; other species not tested.

IHC: 1:20-1:200

IF: 1:10-1:100

Antibody description

Cyclin B2 Rabbit Polyclonal antibody. Positive IF detected in Hela cells. Positive IHC detected in human cervical cancer tissue, human lung cancer tissue, human testis tissue. Positive IP detected in A431 cells. Positive WB detected in A431 cells, K-562 cells. Observed molecular weight by Western-blot: 51 kDa

Preparation

This antibody was obtained by immunization of Cyclin B2 recombinant protein (Accession Number: NM_004701). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, WB, IF, IHC, IP


Dilutions

Recommended Dilution:

Validations




A431 cells were subjected to SDS PAGE followed by western blot with Catalog No:109663(CCNB2 antibody) at dilution of 1:500




Immunofluorescent analysis of Hela cells, using CCNB2 antibody Catalog No:109663 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red).

Cyclin B2 antibody

Catalog Number: 109663



Immunohistochemistry of paraffin-embedded human cervical cancer tissue slide using Catalog No:109663(CCNB2 Antibody) at dilution of 1:200 (under 40x lens)



IP Result of anti-CCNB2 (IP:Catalog No:109663, 4ug; Detection:Catalog No:109663 1:1000) with A431 cells lysate 2000ug.

