

HIST2H2AA4 antibody



Catalog Number: 111363

Product name

HIST2H2AA4 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

HIST2H2AA4 Rabbit Polyclonal antibody. Positive WB detected in L02 cells, Jurkat cells, K-562 cells, mouse liver tissue, rat liver tissue. Positive IHC detected in human liver cancer tissue. Observed molecular weight by Western-blot: 14 kDa

Preparation

This antibody was obtained by immunization of HIST2H2AA4 recombinant protein (Accession Number: BC001629). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IHC, WB


Dilutions

Recommended Dilution:

WB: 1:200-1:2000


IHC: 1:20-1:200

Validations



L02 cells were subjected to SDS PAGE followed by western blot with Catalog

No:111363(HIST2H2AA4 Antibody) at dilution of 1:600



Immunohistochemistry of paraffin-embedded human liver cancer tissue slide using Catalog No:111363(HIST2H2AA4 Antibody) at dilution of 1:50 (under 40x lens)