

ADAM17-Specific antibody



Catalog Number: 107763

Product name

ADAM17-Specific antibody

Specificity

Human; other species not tested.

Antibody description

ADAM17-Specific Rabbit Polyclonal antibody.
Positive WB detected in human placenta tissue.
Positive IHC detected in human pancreas cancer tissue. Observed molecular weight by Western-blot: 29 kDa

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_003183).
Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol
pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200


Validations



human placenta tissue were subjected to SDS PAGE followed by western blot with Catalog No:107763(ADAM17-Specific antibody) at dilution of 1:500



Immunohistochemistry of paraffin-embedded human pancreas cancer tissue slide using Catalog No:107763(ADAM17-Specific Antibody) at dilution of 1:200 (under 10x lens).



Immunohistochemistry of paraffin-embedded human pancreas cancer tissue slide using Catalog No:107763(ADAM17-Specific Antibody) at dilution of 1:200 (under 40x lens).