

# CACNB2 antibody



Catalog Number: 108763

## Product name

CACNB2 antibody

## Immunogen

[Human CACNB2 Recombinant protein \(GST tag & His tag\)](#)

## Specificity

Human, Mouse; other species not tested.

## Antibody description

CACNB2 Rabbit Polyclonal antibody. Positive WB detected in mouse heart tissue. Observed molecular weight by Western-blot: 74 kDa

## Preparation

This antibody was obtained by immunization of CACNB2 recombinant protein (Accession Number: NM\_201596). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

ELISA, WB

## Dilutions

Recommended Dilution:  
WB: 1:200-1:1000

## Validations



mouse heart tissue were subjected to SDS PAGE followed by western blot with Catalog No:108763(CACNB2 Antibody) at dilution of 1:300