

NELFE Polyclonal Antibody



Catalog Number: 164073

Product name

NELFE Polyclonal Antibody

Specificity

Human, Mouse

Antibody description

Polyclonal antibody to NELFE

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-280 of human NELFE (NP_002895.3).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC, IF


Dilutions

WB 1:500 - 1:2000

IHC 1:50 - 1:200


IF 1:50 - 1:200

Validations




Western blot - NELFE Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using NELFE antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 90s.




Immunohistochemistry - NELFE Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human lung cancer using NELFE antibody at dilution of 1:100 (40x lens).



Immunohistochemistry - NELFE Polyclonal Antibody


Immunohistochemistry of paraffin-embedded human colon carcinoma using NELFE antibody at dilution of 1:100 (40x lens).



Immunohistochemistry - NELFE Polyclonal Antibody

Immunohistochemistry of paraffin-embedded


mouse Intestine using NELFE antibody at dilution



of 1:100 (40x lens).

Immunohistochemistry - NELFE Polyclonal Antibody

Immunohistochemistry of paraffin-embedded mouse liver using NELFE antibody at dilution of 1:100 (40x lens).



Immunofluorescence - NELFE Polyclonal Antibody

Immunofluorescence analysis of A-549 cells using NELFE antibody .