

# CKS1B Polyclonal Antibody



Catalog Number: 165073

---

## Product name

CKS1B Polyclonal Antibody

## Specificity

Human, Mouse, Rat

## Antibody description

Polyclonal antibody to CKS1B

## Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-79 of human CKS1B (NP\_001817.1).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

WB

## Dilutions

WB 1:500 - 1:2000

## Validations



Western blot - CKS1B Polyclonal Antibody

Western Blot analysis of Jurkat cell nucleus extracts using Cks1 Polyclonal Antibody.