

# Anti-Hexokinase I (6B8) Mouse antibody



Catalog Number: 168073

## Product name

Anti-Hexokinase I (6B8) Mouse antibody

## Specificity

Human, Mouse, Rat

## Antibody description

Hexokinase I (6B8) Mouse monoclonal antibody

## Preparation

Antigen: Purified recombinant human Hexokinase 1 protein fragments expressed in E.coli

## Formulation

Purified mouse monoclonal antibody in PBS(pH 7.4) containing with 0.02% sodium azide and 50% glycerol.

## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

Mouse IgG1


## Applications

WB

## Dilutions

WB: 1/1000

## Validations



Western blot detection of Hexokinase 1 in MCF7, Jurkat, CHO-K1, C2C12, C6, 3T3 and HeLa cell lysates using Hexokinase 1 mouse mAb (1:1000 diluted). Predicted band size: 102KDa. Observed band size: 102KDa.