

ANKRD27 antibody



Catalog Number: 108073

Product name

ANKRD27 antibody

IP: 1:500-1:5000

IHC: 1:20-1:200

Specificity

Human; other species not tested.

Antibody description

ANKRD27 Rabbit Polyclonal antibody. Positive IHC detected in human breast cancer tissue. Positive WB detected in MCF-7 cells. Positive IP detected in A431 cells, MCF-7 cells. Observed molecular weight by Western-blot: 117-130 kDa

Preparation

This antibody was obtained by immunization of ANKRD27 recombinant protein (Accession Number: XM_006723412). Purification method: Antigen Affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, IHC, IP, WB

Dilutions


Recommended Dilution:

WB: 1:500-1:5000


Validations



MCF-7 cells were subjected to SDS PAGE followed by western blot with Catalog No:108073(ANKRD27 Antibody) at dilution of



IP Result of anti-ANKRD27 (IP:Catalog No:108073, 4ug; Detection:Catalog No:108073 1:1000) with A431 cells lysate 2000ug.



Immunohistochemistry of paraffin-embedded human breast cancer tissue slide using Catalog

ANKRD27 antibody



Catalog Number: 108073

No:108073(ANKRD27 Antibody) at dilution of 1:50 (under 40x lens)