

# SNX18 antibody



Catalog Number: 115473

## Product name

SNX18 antibody

## Immunogen

[Human SNX18 Recombinant protein \(GST tag & His tag\)](#)

## Specificity

Human; other species not tested.

## Antibody description

SNX18 Rabbit Polyclonal antibody. Positive WB detected in HL-60 cells, human placenta tissue, L02 cells, MCF7 cells. Observed molecular weight by Western-blot: 67-70kd

## Preparation

This antibody was obtained by immunization of SNX18 recombinant protein (Accession Number: NM\_001102575). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

## Validations



HL-60 cells were subjected to SDS PAGE followed by western blot with Catalog No:115473(SNX18 antibody) at dilution of 1:1000