

UBAP1 antibody



Catalog Number: 116473

Product name

UBAP1 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

UBAP1 Rabbit Polyclonal antibody. Positive IF detected in Hela cells. Positive IHC detected in human lung cancer tissue, human heart tissue, human lung tissue. Positive IP detected in mouse brain tissue. Positive WB detected in human heart tissue, A431 cells, human lung tissue, human spleen tissue, human testis tissue, mouse brain tissue, mouse lung tissue, mouse testis tissue, mouse thymus tissue, Y79 cells. Observed molecular weight by Western-blot: 55kd

Preparation

This antibody was obtained by immunization of UBAP1 recombinant protein (Accession Number: NM_001171203). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC, IF, IP

Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:200-1:2000


IHC: 1:20-1:200

IF: 1:10-1:100

Validations




human heart tissue were subjected to SDS PAGE followed by western blot with Catalog No:116473(UBAP1 antibody) at dilution of 1:500




Immunohistochemical of paraffin-embedded human lung cancer using Catalog No:116473(UBAP1 antibody) at dilution of 1:50 (under 10x lens)

UBAP1 antibody

Catalog Number: 116473



Immunofluorescent analysis of HeLa cells, using UBAP1 antibody Catalog No:116473 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit



IgG (red).

IP Result of anti-UBAP1 (IP:Catalog No:116473, 4ug; Detection:Catalog No:116473 1:500) with mouse brain tissue lysate 4000ug.