

Anti-TCTP Rabbit antibody



Catalog Number: 166873

Product name

Anti-TCTP Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

TCTP Rabbit polyclonal antibody

Preparation

Antigen: A synthetic peptide of human TCTP

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at -20°C. Stable for 12 months from date of receipt.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, ICC/IF, FC, IP

Dilutions


WB: 1/2000-1/10000

ICC/IF: 1/50


FC: 1/20-1/100

IP: 1/20

Validations



Immunocytochemistry of TCTP(green) in HeLa cells using TCTP Rabbit pAb at dilution 1/50, and



Western blot detection of TCTP in K562,Rat Brain,C6,3T3,HeLa cell lysates using TCTP Rabbit pAb(1:1000 diluted).Predicted band size:20kDa.Observed band size:20kDa.