

EIF1AD antibody



Catalog Number: 110173

Product name

EIF1AD antibody

Immunogen

[Human EIF1AD Recombinant protein \(GST tag & His tag\)](#)

Specificity

human, mouse, rat; other species not tested.

Antibody description

EIF1AD Rabbit Polyclonal antibody. Positive WB detected in K-562 cells, HEK-293 cells, A549 cells, HL-60 cells, MCF7 cells. Observed molecular weight by Western-blot: 27kd

Preparation

This antibody was obtained by immunization of EIF1AD recombinant protein (Accession Number: NM_001242482). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



K-562 cells were subjected to SDS PAGE followed by western blot with Catalog No:110173(EIF1AD antibody) at dilution of 1:500