

Anti-ING1 Rabbit antibody



Catalog Number: 167573

Product name

Rabbit IgG

Anti-ING1 Rabbit antibody

Applications

Specificity

WB, FC

Human, Mouse, Rat

Dilutions

Antibody description

WB: 1/1000

ING1 Rabbit monoclonal antibody

FC: 1/20

Preparation

Antigen: A synthetic peptide of human ING1

Validations

Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl,
40%Glycerol, 0.01% sodium azide and 0.05% BSA


Storage

Store at 4°C short term. Store at -20°C long term.
Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type



Western blot detection of ING1 in HeLa, A549, HL-60, U251, U87-MG cell lysates using ING1 Rabbit mAb(1:1000 diluted). Predicted band size: 47kDa. Observed band size: 47kDa.