

Anti-HIF-1 alpha Rabbit antibody



Catalog Number: 166973

Product name

Anti-HIF-1 alpha Rabbit antibody

IP: 1/30

FC: 1/30

Specificity

Human

Antibody description

HIF-1 alpha Rabbit polyclonal antibody

Preparation

Antigen: A synthesized peptide derived from human HIF-1 alpha

Formulation

Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P, ICC/IF, IP, FC


Dilutions

WB: 1/500-1/2000


IHC: 1/50-1/200

ICC/IF: 1/50-1/200


Validations



Western blot analysis of HIF-1 alpha expression in Ramos cell lysate.



Immunohistochemical analysis of paraffin-embedded human kidney, using HIF-1 alpha



Antibody.

Immunofluorescent analysis of HeLa cells, using HIF-1 alpha Antibody.