

NRBP2 antibody



Catalog Number: 113273

Product name

NRBP2 antibody

Immunogen

[Human NRBP2 Recombinant protein \(GST tag & His tag\)](#)

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

NRBP2 Rabbit Polyclonal antibody. Positive WB detected in mouse kidney tissue, HeLa cells, human placenta tissue, human testis tissue, mouse brain tissue. Positive IP detected in mouse kidney tissue. Observed molecular weight by Western-blot: 33 kDa, 58 kDa

Preparation

This antibody was obtained by immunization of NRBP2 recombinant protein (Accession Number: BC113873). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:200-1:2000

Validations



mouse kidney tissue were subjected to SDS PAGE followed by western blot with Catalog No:113273(NRBP2 antibody) at dilution of 1:500



IP Result of anti-NRBP2 (IP:Catalog No:113273, 4ug; Detection:Catalog No:113273 1:500) with mouse kidney tissue lysate 4000ug.