

Frizzled 7 antibody



Catalog Number: 110773

Product name

Frizzled 7 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

Frizzled 7 Rabbit Polyclonal antibody. Positive IP detected in mouse skeletal muscle tissue. Positive WB detected in mouse skeletal muscle tissue, mouse heart tissue, mouse kidney tissue. Observed molecular weight by Western-blot: 64 kDa

Preparation

This antibody was obtained by immunization of Frizzled 7 recombinant protein (Accession Number: NM_003507). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:200-1:2000

Validations



mouse skeletal muscle tissue were subjected to SDS PAGE followed by western blot with Catalog No:110773(FZD7 antibody) at dilution of 1:600



IP Result of anti-FZD7 (IP:Catalog No:110773, 4ug; Detection:Catalog No:110773 1:600) with mouse skeletal muscle tissue lysate 2200ug.