

Anti-NELFe Rabbit antibody



Catalog Number: 167083

Product name

WB: 1/1000

Anti-NELFe Rabbit antibody

IHC: 1/20

Specificity

ICC/IF: 1/20

Human, Mouse, Rat

IP: 1/20

Antibody description

NELFe Rabbit polyclonal antibody

Preparation

Antigen: Recombinant protein of human NELFe

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at -20°C. Stable for 12 months from date of receipt.

Clonality

Polyclonal

Ig Type


Rabbit IgG

Applications


WB, IHC-P, ICC/IF, IP

Dilutions

Validations



Immunohistochemistry of NELFe in paraffin-embedded Human breast cancer tissue using



NELFe Rabbit pAb at dilution 1/20

Western blot detection of NELFe in C6 cell lysates using NELFe Rabbit pAb(1:1000 diluted). Predicted band size:43kDa.Observed band size:43kDa.