

TBRG1 antibody



Catalog Number: 115883

Product name

TBRG1 antibody

Immunogen

[Human TBRG1 Recombinant protein \(His tag\)](#)

Specificity

human, mouse, rat; other species not tested.

Antibody description

TBRG1 Rabbit Polyclonal antibody. Positive WB detected in BxPC-3 cells. Observed molecular weight by Western-blot: 45kd

Preparation

This antibody was obtained by immunization of TBRG1 recombinant protein (Accession Number: BC032312). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, WB

Dilutions

Recommended Dilution:
WB: 1:200-1:2000

Validations



BxPC-3 cells were subjected to SDS PAGE followed by western blot with Catalog No:115883(TBRG1 antibody) at dilution of 1:600