

FUT8 antibody



Catalog Number: 107283

Product name

FUT8 antibody

IHC: 1:20-1:200


Specificity

Human; other species not tested.

Validations

Antibody description

FUT8 Mouse Monoclonal antibody. Positive IHC detected in human lung cancer tissue, human breast cancer tissue. Positive WB detected in HeLa cells, PC-3 cells. Observed molecular weight by Western-blot: 66 kDa




Preparation

This antibody was obtained by immunization of FUT8 recombinant protein (Accession Number: BC025385). Purification method: Protein A purified.

HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:107283(FUT8 antibody) at dilution of 1:1000

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.



Storage

Store at -20°C. DO NOT ALIQUOT


Clonality

Monoclonal

Immunohistochemistry of paraffin-embedded human breast cancer slide using Catalog No:107283(FUT8 Antibody) at dilution of 1:50 under 40x lens

Ig Type

Mouse IgG2a



Applications

ELISA, WB, IHC

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Immunohistochemistry of paraffin-embedded human lung cancer tissue slide using Catalog No:107283(FUT8 Antibody) at dilution of 1:50

FUT8 antibody



Catalog Number: 107283

under 40x lens