

DNAJB2 antibody



Catalog Number: 110093

Product name

DNAJB2 antibody

Immunogen

[Human DNAJB2 Recombinant protein \(GST tag & His tag\)](#)

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

DNAJB2 Rabbit Polyclonal antibody. Positive WB detected in K-562 cells, A431 cells, HeLa cells, MCF7 cells, mouse brain tissue, mouse heart tissue, mouse skeletal muscle tissue. Positive IP detected in mouse brain tissue. Observed molecular weight by Western-blot: 36-40 kDa

Preparation

This antibody was obtained by immunization of DNAJB2 recombinant protein (Accession Number: NM_006736). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:500-1:5000

Validations



K-562 cells were subjected to SDS PAGE followed by western blot with Catalog No:110093(DNAJB2 antibody) at dilution of 1:300



IP Result of anti-DNAJB2 (IP:Catalog No:110093, 4ug; Detection:Catalog No:110093 1:1000) with mouse brain tissue lysate 4000ug.