

Anti-CHK2/CHEK2 antibody



Catalog Number: 104093

Product name

Anti-CHK2/CHEK2 antibody

Immunogen

[Mouse CHK2/CHEK2 \(His & GST Tag\) recombinant protein](#)

Specificity

Mouse CHK2 / CHEK2

Antibody description

Rabbit polyclonal to CHK2/CHEK2

Preparation

Produced in rabbits immunized with purified, recombinant Mouse CHK2 / CHEK2 (rM CHK2 / CHEK2; Q9Z265; Mey1-Leu 546). CHK2 / CHEK2 specific IgG was purified by Mouse CHK2 / CHEK2 affinity chromatography.

Formulation

0.2 µm filtered solution in PBS with 5% trehalose

Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P, IP


Dilutions

WB: 10-20 µg/mL


IHC-P: 0.1-2 µg/mL

IP: 4-6 µg/mg of lysate

Validations



CHK2 / CHEK2 Antibody, Rabbit PAb, Antigen Affinity Purified, Western blot




CHK2 / CHEK2 Antibody, Rabbit PAb, Antigen Affinity Purified, Immunoprecipitation

Anti-CHK2/CHEK2 antibody



Catalog Number: 104093



CHK2 / CHEK2 Antibody, Rabbit PAb, Antigen Affinity Purified, Immunohistochemistry

Immunochemical staining of mouse CHEK2 in mouse spleen with rabbit polyclonal antibody (1 $\mu\text{g}/\text{mL}$, formalin-fixed paraffin embedded sections).