

# TDG Polyclonal Antibody



Catalog Number: 161593

---

## Product name

TDG Polyclonal Antibody

## Specificity

Human, Mouse, Rat

## Antibody description

Polyclonal antibody to TDG

## Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 141-410 of human TDG (NP\_003202.3).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


WB, IF

## Dilutions

WB 1:500 - 1:2000

IF 1:50 - 1:200

## Validations



Immunofluorescence - TDG Polyclonal Antibody

Immunofluorescence analysis of A549 cells using TDG antibody . Blue: DAPI for nuclear staining.