

ICK antibody



Catalog Number: 111593

Product name

ICK antibody

Polyclonal

Ig Type

Rabbit IgG

Specificity

human, mouse, rat; other species not tested.

Applications

ELISA, WB

Antibody description

ICK Rabbit Polyclonal antibody. Positive WB detected in PC-3 cells. Observed molecular weight by Western-blot: 71kd

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Preparation

This antibody was obtained by immunization of ICK recombinant protein (Accession Number: NM_173041). Purification method: Antigen affinity purified.

Validations



Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

PC-3 cells were subjected to SDS PAGE followed by western blot with Catalog No:111593(ICK antibody) at dilution of 1:300

Clonality