

Anti-Phospho-Nrf2 (Ser40) Rabbit antibody



Catalog Number: 168593

Product name

Rabbit IgG

Anti-Phospho-Nrf2 (Ser40) Rabbit antibody

Applications

Specificity

WB, IHC-P, ICC/IF, FC

Human

Dilutions

Antibody description

1:2000

Phospho-Nrf2 (Ser40) Rabbit monoclonal antibody

1:50

Preparation

1:50

Antigen: A synthetic phosphopeptide corresponding to residues surrounding Ser40 of human Nrf2

1:20

Validations

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.


Storage

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type



Western blot detection of Phospho-Nrf2 (Ser40) in K562, HeLa cell lysates using Phospho-Nrf2 (Ser40) Rabbit mAb(1:1000 diluted). Predicted band size:68kDa. Observed band size:100kDa.