

Anti-GNAQ Rabbit antibody



Catalog Number: 167293

Product name

Rabbit IgG

Anti-GNAQ Rabbit antibody

Applications

Specificity

WB, IHC-P, FC, IP

Human, Mouse, Rat

Dilutions

Antibody description

WB: 1/1000

GNAQ Rabbit polyclonal antibody

IHC: 1/200

Preparation

FC: 1/20

Antigen: Recombinant protein of human GNAQ

IP: 1/20

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage


Store at -20°C. Stable for 12 months from date of receipt.

Clonality

Polyclonal

Ig Type

Validations



Western blot detection of GNAQ in K562,Rat Brain,C6,Hela cell lysates using GNAQ Rabbit pAb(1:1000 diluted).Predicted band size:42kDa.Observed band size:42kDa.