

Anti-ELOVL5 Rabbit antibody



Catalog Number: 168693

Product name

WB, ICC/IF, FC

Anti-ELOVL5 Rabbit antibody

Dilutions

Specificity

1:1000

Human, Mouse, Rat

1:20

Antibody description

1:20

ELOVL5 Rabbit monoclonal antibody

Preparation

Antigen: A synthetic peptide of human ELOVL5

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze / thaw cycle.

Clonality


Monoclonal

Ig Type


Rabbit IgG

Applications

Validations



Immunocytochemistry of ELOVL5(green) in HeLa cells using ELOVL5 Rabbit mAb at dilution 1/200, and DAPI(blue)



and DAPI(blue)

Western blot detection of ELOVL5 in Jurkat, HeLa cell lysates using ELOVL5 Rabbit mAb(1:1000 diluted). Predicted band size: 35kDa. Observed band size: 35kDa.