

HEATR4 antibody



Catalog Number: 111393

Product name

HEATR4 antibody

Immunogen

[Human HEATR4 Recombinant protein \(His tag\)](#)

Specificity

Human; other species not tested.

Antibody description

HEATR4 Rabbit Polyclonal antibody. Positive IHC detected in human testis tissue, human breast cancer tissue.

Preparation

This antibody was obtained by immunization of HEATR4 recombinant protein (Accession Number: BC047590). Purification method: Antigen Affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, IHC

Dilutions


Recommended Dilution:

IHC: 1:20-1:200

Validations



Immunohistochemistry of paraffin-embedded human testis slide using Catalog No:111393(HEATR4 Antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemistry of paraffin-embedded human testis slide using Catalog No:111393(HEATR4 Antibody) at dilution of 1:50 (under 40x lens)