

Anti-Heme Oxygenase 1 Rabbit antibody



Catalog Number: 167393

Product name

Anti-Heme Oxygenase 1 Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

Heme Oxygenase 1 Rabbit polyclonal antibody

Preparation

Antigen: Recombinant protein of mouse Heme Oxygenase 1

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at -20°C. Stable for 12 months from date of receipt.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P, ICC/IF, FC, IP

Dilutions

WB: 1/1000


IHC: 1/50

ICC/IF: 1/50


FC: 1/20

IP: 1/20


Validations



Immunohistochemistry of Heme Oxygenase 1 in paraffin-embedded Human lung cancer tissue using Heme Oxygenase 1 Rabbit pAb at dilution



Western blot detection of Heme Oxygenase 1 in U251MG,K562,A549 cell lysates using Heme Oxygenase 1 Rabbit pAb(1:500 diluted).Predicted band size:33kDa.Observed band size:33kDa.



Western blot detection of Heme Oxygenase 1 in Mouse spleen lysates using Heme Oxygenase 1

Anti-Heme Oxygenase 1 Rabbit antibody



Catalog Number: 167393

Rabbit pAb(1:1000 diluted).Predicted band

size:33KDa.Observed band size:33KDa.