

CLIP-115/CYLN2 antibody



Catalog Number: 109393

Product name

CLIP-115/CYLN2 antibody

Immunogen

[Human CLIP2 Recombinant protein \(GST tag & His tag\)](#)

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

CLIP-115/CYLN2 Rabbit Polyclonal antibody. Positive WB detected in mouse brain tissue, rat brain tissue. Observed molecular weight by Western-blot: 110-120 kDa

Preparation

This antibody was obtained by immunization of CLIP-115/CYLN2 recombinant protein (Accession Number: NM_032421). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:1000-1:10000

Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:109393(CLIP2 Antibody) at dilution of 1:3000