

ING3 antibody



Catalog Number: 111793

Product name

ING3 antibody

Immunogen

[Human ING3 Recombinant protein \(GST tag & His tag\)](#)

Specificity

Human, Mouse; other species not tested.

Antibody description

ING3 Rabbit Polyclonal antibody. Positive IHC detected in human breast cancer tissue. Positive WB detected in mouse uterus tissue, HepG2 cells. Observed molecular weight by Western-blot: 47 kDa

Preparation

This antibody was obtained by immunization of ING3 recombinant protein (Accession Number: NM_198267). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IHC, WB


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200

Validations



mouse uterus tissue were subjected to SDS PAGE followed by western blot with Catalog No:111793(ING3 antibody) at dilution of 1:800



Immunohistochemical of paraffin-embedded human breast cancer using Catalog No:111793(ING3 antibody) at dilution of 1:100 (under 10x lens)