

VPS52 antibody



Catalog Number: 116793

Product name

IP: 1:500-1:5000

VPS52 antibody

IHC: 1:20-1:200

Specificity

IF: 1:10-1:100

Human, Mouse, Rat; other species not tested.

Validations

Antibody description

VPS52 Rabbit Polyclonal antibody. Positive IF detected in HepG2 cells. Positive IHC detected in human gliomas tissue. Positive IP detected in mouse skin tissue. Positive WB detected in HepG2 cells, COLO 320 cells, mouse skin tissue. Observed molecular weight by Western-blot: 82kd

Preparation

This antibody was obtained by immunization of VPS52 recombinant protein (Accession Number: NM_022553). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, WB, IHC, IF, IP


Dilutions

Recommended Dilution:


WB: 1:200-1:2000



HepG2 cells were subjected to SDS PAGE followed by western blot with Catalog No:116793(VPS52 antibody) at dilution of 1:1000



Immunohistochemical of paraffin-embedded human gliomas using Catalog No:116793(VPS52 antibody) at dilution of 1:50 (under 10x lens)




VPS52 antibody



Catalog Number: 116793

IP Result of anti-VPS52 (IP:Catalog No:116793, 3ug; Detection:Catalog No:116793 1:1000) with mouse skin tissue lysate 3200ug.



Immunofluorescent analysis of HepG2 cells using Catalog No:116793(VPS52 Antibody) at dilution of 1:25 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)