

GNL1 antibody



Catalog Number: 111004

Product name

GNL1 antibody

IHC: 1:20-1:200

IF: 1:10-1:100

Specificity

human; other species not tested.

Antibody description

GNL1 Rabbit Polyclonal antibody. Positive IHC detected in human heart. Positive IF detected in Hela. Positive WB detected in MCF7 cells, HeLa cells, Jurkat cells, K-562 cells, Raji cells, U-937 cells. Observed molecular weight by Western-blot: 69-75kd

Preparation

This antibody was obtained by immunization of GNL1 recombinant protein (Accession Number: XM_005249015). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB, IHC, IF

Dilutions


Recommended Dilution:

WB: 1:500-1:5000


Validations



MCF7 cells were subjected to SDS PAGE followed by western blot with Catalog No:111004(GNL1 antibody) at dilution of 1:300



Immunohistochemical of paraffin-embedded human heart using Catalog No:111004(GNL1 antibody) at dilution of 1:100 (under 10x lens)




Immunohistochemical of paraffin-embedded human heart using Catalog No:111004(GNL1

GNL1 antibody



Catalog Number: 111004

antibody) at dilution of 1:100 (under 40x lens)



Immunofluorescent analysis of HeLa cells, using GNL1 antibody Catalog No:111004 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red).