

MYT1 antibody



Catalog Number: 113004

Product name

MYT1 antibody

Immunogen

[Human MYT1 Recombinant protein \(GST tag & His tag\)](#)

Specificity

Human; other species not tested.

Antibody description

MYT1 Rabbit Polyclonal antibody. Positive WB detected in SH-SY5Y cells. Observed molecular weight by Western-blot: 122 kDa

Preparation

This antibody was obtained by immunization of MYT1 recombinant protein (Accession Number: BC062313). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, WB

Dilutions

Recommended Dilution:
WB: 1:200-1:1000

Validations



SH-SY5Y cells were subjected to SDS PAGE followed by western blot with Catalog No:113004(MYT1 Antibody) at dilution of 1:300