

Anti-UBE2A antibody



Catalog Number: 176404

Product name

Anti-UBE2A antibody

Specificity

Human, Mouse, Rat

Antibody description

Rabbit polyclonal antibody to UBE2A

Preparation

This antigen of this antibody was synthesized peptide derived from the c-terminal region of human ube2a/b.

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P

Dilutions

WB: 1:500-1:2000

IHC-P: 1:100-300

Validations




Fig1.; Western Blot analysis of extracts from Jurkat cells, using UBE2A/B Polyclonal Antibody.. Secondary antibody (catalog#: HA1001) was diluted at 1:20000




Fig2.; Immunohistochemical analysis of paraffin-embedded rat-brain, antibody was diluted at 1:100