

INO80 antibody



Catalog Number: 111804

Product name

INO80 antibody

Immunogen

[Human INO80 Recombinant protein \(GST tag & His tag\)](#)

Specificity

Human; other species not tested.

Antibody description

INO80 Rabbit Polyclonal antibody. Positive IP detected in HeLa cells. Positive WB detected in HeLa cells, HepG2 cells. Observed molecular weight by Western-blot: 177,150kDa

Preparation

This antibody was obtained by immunization of INO80 recombinant protein (Accession Number: NM_017553). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:2000

Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:111804(INO80 Antibody) at dilution of 1:600



IP Result of anti-INO80 (IP:Catalog No:111804, 4ug; Detection:Catalog No:111804 1:500) with HeLa cells lysate 2000ug.