

# RPGRIP1 antibody



Catalog Number: 114804

## Product name

RPGRIP1 antibody

## Immunogen

[Human RPGRIP1 Recombinant protein \(GST tag & His tag\)](#)

## Specificity

Human; other species not tested.

## Antibody description

RPGRIP1 Rabbit Polyclonal antibody. Positive WB detected in HEK-293 cells. Observed molecular weight by Western-blot: 70 kDa, 140-150 kDa

## Preparation

This antibody was obtained by immunization of RPGRIP1 recombinant protein (Accession Number: NM\_001377523). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

ELISA, WB

## Dilutions

Recommended Dilution:  
WB: 1:500-1:5000

## Validations



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:114804(RPGRIP1 antibody) at dilution of 1:400