

Anti-PSMD14 Rabbit antibody



Catalog Number: 166904

Product name

Anti-PSMD14 Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

PSMD14 Rabbit polyclonal antibody

Preparation

Antigen: A synthetic peptide of human PSMD14

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at -20°C. Stable for 12 months from date of receipt.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P, ICC/IF


Dilutions

WB: 1/1000-1/5000


IHC: 1/50-1/200

ICC/IF: 1/50-1/200

Validations



Immunocytochemistry of PSMD14(green) in HeLa cells using PSMD14 Rabbit pAb at dilution 1/50, and DAPI(blue)



and DAPI(blue)

Western blot detection of PSMD14 in HeLa, COS7, PC12 cell lysates using PSMD14 Rabbit pAb(1:1000 diluted). Predicted band size: 35KDa. Observed band size: 35KDa.