

FAS/CD95 antibody



Catalog Number: 107204

Product name

FAS/CD95 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

FAS/CD95 Mouse Monoclonal antibody. Positive FC detected in Jurkat cells. Positive WB detected in MCF7 cells. Observed molecular weight by Western-blot: 38 kDa

Preparation

This antibody was obtained by immunization of FAS/CD95 recombinant protein (Accession Number: NM_000043). Purification method: Protein G purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type

Mouse IgG1

Applications


ELISA, WB, FC

Dilutions


Recommended Dilution:

WB: 1:200-1:1000

Validations



MCF7 cells were subjected to SDS PAGE followed by western blot with Catalog No:107204(FAS antibody) at dilution of 1:500



1×10^6 Jurkat cells were stained with 0.2ug FAS antibody (Catalog No:107204, red) and control antibody (blue). Fixed with 90% MeOH blocked with 3% BSA (30 min). Alexa Fluor 488-conjugated AffiniPure Goat Anti-Mouse IgG(H+L) with dilution 1:1000.