

GC Polyclonal Antibody



Catalog Number: 164604

Product name

GC Polyclonal Antibody

Specificity

Human, Mouse, Rat

Antibody description

Polyclonal antibody to GC

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 324-493 of human GC (NP_001191236.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IF

Dilutions

WB 1:500 - 1:2000


IF 1:10 - 1:100

Validations



Western blot - GC Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using GC antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure



Immunofluorescence - GC Polyclonal Antibody

Immunofluorescence analysis of U2OS cells using GC antibody. Blue: DAPI for nuclear staining.