

RBM19 antibody



Catalog Number: 114604

Product name

RBM19 antibody

IHC: 1:20-1:200

IF: 1:10-1:100

Specificity

Human, Mouse; other species not tested.

Antibody description

RBM19 Rabbit Polyclonal antibody. Positive WB detected in mouse colon tissue, COLO 320 cells, mouse small intestine tissue. Positive IF detected in MCF-7 cells. Positive IHC detected in human colon tissue. Observed molecular weight by Western-blot: 120kd

Preparation

This antibody was obtained by immunization of RBM19 recombinant protein (Accession Number: NM_001146698). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB, IHC, IF

Dilutions


Recommended Dilution:

WB: 1:500-1:5000


Validations



mouse colon tissue were subjected to SDS PAGE followed by western blot with Catalog No:114604(RBM19 antibody) at dilution of 1:300



Immunohistochemical of paraffin-embedded human colon using Catalog No:114604(RBM19 antibody) at dilution of 1:50 (under 10x lens)




Immunohistochemical of paraffin-embedded human colon using Catalog No:114604(RBM19

RBM19 antibody



Catalog Number: 114604

antibody) at dilution of 1:50 (under 40x lens)



Immunofluorescent analysis of MCF-7 cells, using RBM19 antibody Catalog No:114604 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red).