

Ch-TOG antibody



Catalog Number: 109304

Product name

Ch-TOG antibody

Immunogen

[Human CKAP5 Recombinant protein \(GST tag & His tag\)](#)

Specificity

Human; other species not tested.

Antibody description

Ch-TOG Rabbit Polyclonal antibody. Positive WB detected in NCCIT cell. Observed molecular weight by Western-blot: 220 kDa

Preparation

This antibody was obtained by immunization of Ch-TOG recombinant protein (Accession Number: NM_001008938). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, WB

Dilutions

Recommended Dilution:
WB: 1:500-1:5000

Validations



NCCIT cell were subjected to SDS PAGE followed by western blot with Catalog No:109304(CKAP5 Antibody) at dilution of 1:1000