

# Amylase, alpha antibody



Catalog Number: 108014

## Product name

Amylase, alpha antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

Amylase, alpha Rabbit Polyclonal antibody.  
Positive IHC detected in human pancreas tissue.  
Positive WB detected in human liver tissue, mouse pancreas tissue. Observed molecular weight by Western-blot: 58 kDa

## Preparation

This antibody was obtained by immunization of Amylase, alpha recombinant protein (Accession Number: NM\_000699). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200


## Validations



human liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:108014(AMY2A antibody) at dilution of 1:500



Immunohistochemical of paraffin-embedded human pancreas using Catalog No:108014(AMY2A antibody) at dilution of 1:100 (under 10x lens)



Immunohistochemical of paraffin-embedded human pancreas using Catalog No:108014(AMY2A antibody) at dilution of 1:100 (under 40x lens)