

TRPS1 antibody



Catalog Number: 116414

Product name

TRPS1 antibody

IP: N/AIHC: 1:20-1:200

IF: 1:20-1:200

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

TRPS1 Rabbit Polyclonal antibody. Positive IP detected in HeLa cells. Positive WB detected in mouse brain tissue, mouse kidney tissue, rat brain tissue. Positive IF detected in HepG2 cells. Positive IHC detected in human prostate cancer tissue, human kidney tissue. Observed molecular weight by Western-blot: 135-142kd

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_014112).
Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol
pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB, IF, IHC, IP

Dilutions


Recommended Dilution:

WB: 1:200-1:2000

Validations




mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:116414(TRPS1 antibody) at dilution of 1:100



Immunoprecipitation
2nd Antibody: Light and Heavy chain
IP-TRPS1: IP of HeLa cell lysate(500ug)
MW of Target band: 142-175kd


IP result of anti-TRPS1(Catalog No:116414 for IP and Detection).




Immunohistochemistry of paraffin-embedded human prostate cancer slide using Catalog No:116414(TRPS1 Antibody) at dilution of 1:50

TRPS1 antibody

Catalog Number: 116414



Immunohistochemistry of paraffin-embedded human prostate cancer slide using Catalog No:116414(TRPS1 Antibody) at dilution of 1:50



Immunofluorescent analysis of HepG2 cells, using TRPS1 antibody Catalog No:116414 at 1:50 dilution and Rhodamine-labeled goat anti-rabbit IgG (red).