

# PHF5A antibody



Catalog Number: 113814

## Product name

PHF5A antibody

## Immunogen

[Human INI-1 Recombinant protein \(GST tag & His tag\)](#)

## Specificity

Human, Mouse; other species not tested.

## Antibody description

PHF5A Rabbit Polyclonal antibody. Positive IP detected in HeLa cells. Positive WB detected in HeLa cells, HEK-293 cells. Observed molecular weight by Western-blot: 12-14 kDa

## Preparation

This antibody was obtained by immunization of PHF5A recombinant protein (Accession Number: NM\_032758). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:2000

## Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:113814(PHF5A antibody) at dilution of 1:500



IP Result of anti-INI-1 (IP:Catalog No:113814, 3ug; Detection:Catalog No:113814 1:500) with HeLa cells lysate 4000ug.