

Anti-PFAS antibody



Catalog Number: 175114

Product name

Anti-PFAS antibody

Specificity

Human

Antibody description

Rabbit polyclonal antibody to PFAS

Preparation

This antigen of this antibody was recombinant protein within human purl aa 1173-1338.

Formulation

Liquid, 1*PBS (pH7.4), 0.2% BSA, 50% Glycerol.
Preservative: 0.05% Sodium Azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P

Dilutions

WB: 1:500

IHC-P: 1:50-1:200

Validations





Fig1: Western blot analysis of PFAS on Hela cell lysate using anti-PFAS antibody at 1/1,000



dilution.

Fig2: Immunohistochemical analysis of paraffin-embedded human tonsil tissue using anti-PFAS antibody. Counter stained with hematoxylin.




Fig3: Immunohistochemical analysis of paraffin-embedded human spleen tissue using anti-PFAS antibody. Counter stained with hematoxylin.




Fig4: Immunohistochemical analysis of paraffin-embedded human breast tissue using anti-PFAS antibody. Counter stained with hematoxylin.




Fig5: Immunohistochemical analysis of paraffin-embedded human kidney tissue using anti-PFAS antibody. Counter stained with hematoxylin.