

# TCOF1 antibody



Catalog Number: 115914

## Product name

TCOF1 antibody

Recommended Dilution:

WB: 1:500-1:5000

## Immunogen

[Human TCOF1 Recombinant protein \(GST tag & His tag\)](#)

IHC: 1:20-1:200

IF: 1:20-1:200

## Specificity

Human; other species not tested.

## Antibody description

TCOF1 Rabbit Polyclonal antibody. Positive IHC detected in human lymphoma tissue. Positive IF detected in MCF-7 cells. Positive WB detected in HeLa cells, Jurkat cells. Observed molecular weight by Western-blot: 220 kDa

## Preparation

This antibody was obtained by immunization of TCOF1 recombinant protein (Accession Number: NM\_001008657). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type


Rabbit IgG

## Applications


ELISA, WB, IHC, IF

## Dilutions

## Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:115914(TCOF1 antibody) at dilution of 1:1000




Immunohistochemical of paraffin-embedded human lymphoma using Catalog No:115914(TCOF1 antibody) at dilution of 1:50 (under 10x lens)

# TCOF1 antibody



Catalog Number: 115914

---



Immunofluorescent analysis of MCF-7 cells, using TCOF1 antibody Catalog No:115914 at 1:50 dilution and FITC-labeled donkey anti-rabbit IgG(green). Blue pseudocolor = DAPI (fluorescent DNA dye).