

# Anti-PBR/TSPO Rabbit antibody



Catalog Number: 168914

## Product name

Anti-PBR/TSPO Rabbit antibody

## Applications

WB, IHC-P, ICC/IF, FC, IP

## Specificity

Human, Mouse

## Dilutions

1:2000

## Antibody description

PBR/TSPO Rabbit monoclonal antibody

1:50

1:50

## Preparation

Antigen: A synthetic peptide of human PBR

1:20

1:20

## Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

## Validations

## Storage


Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

Rabbit IgG



Western blot detection of PBR/TSPO in 3T3, HeLa cell lysates using PBR/TSPO Rabbit mAb(1:1000 diluted). Predicted band size: 19kDa. Observed band size: 19kDa.